



# TG1 Electrocompetent Cell Preparation For Phage Display Library Construction

## Materials

- TG1 *E. coli* strain
  - SOB medium (see preparation below)
  - LB agar plates (see preparation below)
  - Sterile 10% glycerol
  - Sterile ddH<sub>2</sub>O
  - Sterile centrifuge bottles and microfuge tubes
  - Electroporation cuvettes (0.1 cm gap)
  - Electroporator
  - Shaker incubator (37°C)
  - Ice bucket
- 

## Protocol Steps

### Day 1: Starter Culture

1. Streak TG1 cells on LB agar plate.
2. Incubate overnight at 37°C.
3. Inoculate 10 mL SOB medium with a single colony.
4. Grow overnight at 37°C, 200–250 rpm.

### Day 2: Bulk Culture

1. Inoculate 1 L SOB with 10 mL overnight culture (1:100 dilutions).
2. Grow at 37°C, 200 rpm until OD<sub>600</sub> = **0.4–0.5**.
3. Chill culture on ice for 15–30 minutes.

### Cell Harvesting

1. Centrifuge at 4°C, 3000 × g for 10 min.
2. Discard supernatant, resuspend pellet in ice-cold 10% glycerol.
3. Repeat wash **3 times** with 10% glycerol.



4. Final resuspension: ~1 mL per 100 mL original culture.

### Aliquot and Storage

1. Aliquot 50  $\mu$ L into sterile microfuge tubes.
  2. Snap-freeze in liquid nitrogen.
  3. Store at  $-80^{\circ}\text{C}$ .
- 

### Electroporation Conditions

- DNA: 10 pg control plasmid (e.g., pUC19)
  - Voltage: 1.8 kV
  - Capacitance: 25  $\mu$ F
  - Resistance: 200  $\Omega$
  - Time constant: ~5 ms
  - Recovery: Add 950  $\mu$ L SOC, incubate 1 hr at  $37^{\circ}\text{C}$
- 

### Expected Efficiency

- $1-4 \times 10^{10}$  cfu/ $\mu$ g DNA with optimal conditions.
- 

## Media Preparation

### LB Agar Plates

#### Ingredients:

- Tryptone: 10 g
- Yeast extract: 5 g
- NaCl: 10 g
- Agar: 15 g
- Distilled water: 1 L

#### Preparation:



1. Combine tryptone, yeast extract, NaCl, and agar in 1 L distilled water.
2. Heat with stirring until agar dissolves completely.
3. Autoclave at 121°C for 15 minutes.
4. Cool to ~50°C, pour into sterile petri dishes.
5. Allow to solidify at room temperature.
6. Store plates at 4°C, wrapped to prevent drying.

## SOB Medium

### Ingredients:

- Tryptone: 20 g
- Yeast extract: 5 g
- NaCl: 0.5 g
- KCl: 0.186 g
- MgCl<sub>2</sub>·2H<sub>2</sub>O: 10 mL of 1 M stock (final 10 mM)
- MgSO<sub>4</sub>·2H<sub>2</sub>O: 10 mL of 1 M stock (final 10 mM)
- Distilled water: 1 L

### Preparation:

1. Dissolve tryptone, yeast extract, NaCl, and KCl in ~900 mL distilled water.
2. Adjust volume to 1 L with distilled water.
3. Autoclave at 121°C for 15 minutes.
4. Cool to room temperature.
5. Add sterile MgCl<sub>2</sub> and MgSO<sub>4</sub> stock solutions.
6. Mix well before use.