



Anti-fluorescein (FITC) Mouse Monoclonal Antibody

SKU: EAB-004

Recombinant mouse monoclonal antibody (Clone ID: 4M5.3), expressed in Chinese Hamster Ovary cells (CHO), is capable of binding to fluorescein.

species reactivity	independent
recombinant	expressed in Chinese Hamster Ovary cells (CHO).
applications	WB, Flow Cyt, ELISA
antibody form	affinity purified immunoglobulin
immunogen	hapten fluorescein
clone	4M5.3
purity	>95% (SDS-PAGE)
form	0.015 M PBS, 0.05% NaN ₃ , pH7.2
concentration	~ 2 mg/ml
isotype	mouse IgG2a, k
Fc-engineered	n/a; wild type

- Store at -20°C. Recombinant monoclonal antibodies are guaranteed stable for 12 months when properly stored.

References:

1. Kranz, D. M. & Voss, E. W., Jr (1981). Partial elucidation of an anti-hapten repertoire in BALB/c mice: comparative characterization of several mono-clonal anti-fluorescyl antibodies. *Mol. Immunol.* 18, 889–898.
2. Kranz, D. M., Herron, J. N. & Voss, E. W., Jr (1982). Mechanisms of ligand binding by monoclonal anti-fluorescyl antibodies. *J. Biol. Chem.* 257, 6987–6995.
3. Midelfort, K.S., Hernandez, H.H., Lippow, S.M., Tidor, B., Drennan, C.L. and Wittrup, K.D. (2004). Substantial Energetic Improvement with Minimal Structural Perturbation in a High Affinity Mutant Antibody. *J. Mol. Biol.* 343, 685–701.